


AlcoHAWK® Beacon Quick Guide



Be sure to read the entire owner's manual and follow procedures carefully before using this device.

Wait 20 minutes after the subject smokes, eats or drinks before testing. Failure to wait 20 minutes can give inaccurate BAC readings and harm the sensor.

Operating Instructions

- Step 1)** Place a fitted cover over the mouthpiece each time you test another individual.
- Step 2)** If it has been more than 24 hours since the device was last used, **prime the sensor** by blowing into the mouthpiece for 4 to 5 seconds **before** turning on the device.

Step 3) **PRESS** the Power Button and **HOLD** for **one (1) second**. A single beep will sound and the device will display the room temperature.


Note: The device will shut off after 30 seconds of non-use.

Step 4) **PRESS** the Left Button and the device will count down from 99 as it prepares the sensor for testing.

Note: The device will restart a count down from 10 if the sensor is not ready within the first countdown.

Step 5) **WAIT** until the unit displays "rdY" and you hear a **SINGLE BEEP**. Take a **DEEP BREATH AND BLOW STEADILY** for at least 5 seconds. The unit will display "bL" and count from 1 to 5 while you blow. A **SINGLE BEEP** sounds to signal completion of the testing. Your estimated BAC will be displayed in the LCD and the backlight indicates the alert level.

Note: The device will shut off after 10 seconds of non-use.



DO NOT DRINK ALCOHOL AND DRIVE

SEE OWNER'S MANUAL FOR HELP UNDERSTANDING THE RESULTS, PRECAUTIONS AND RE-CALIBRATION INSTRUCTIONS FOR RE-CALIBRATION THAT MUST BE ROUTINELY PERFORMED.

THIS DEVICE IS INTENDED TO MEASURE ALCOHOL IN THE HUMAN BREATH.

Calibration Log

Date First Used	Recalibration Date	Notes