

D I G I T A L

Alcohol Detector

AlcoHAWK[®] Slim

Quick Results

**Easy
To Use**

**FDA 510(k)
Pre-Marketing Notification**

**Meets DOT/NHSTA
Standards for Screeners**

Sleek Design

**2 Digit
Display**

**Multi-User
Mouthpiece Design**



- The AlcoHAWK[®] Slim offers a sleek and portable digital breath alcohol screener at an affordable price. The folding mouthpiece design with replaceable mouthpieces allows for testing of multiple individuals.
- The AlcoHAWK[®] Slim contains a high grade semi-conductor sensor for fast testing for the presence of alcohol.
- The sleek design makes the AlcoHAWK[®] Slim an excellent choice for personal use.

AlcoHAWK[®] Slim

Alcohol Detector *AlcoHAWK Slim*

The AlcoHAWK® Slim is a breath alcohol screener that quickly screens for the presence of alcohol. The AlcoHAWK® Slim is compact, lightweight and offers an innovative folding mouthpiece design. Simply blow into the unit and within seconds the LED displays a digital readout. The AlcoHAWK® Slim includes Sensor Recognition Technology™ to monitor certain types of abnormal sensor conditions. Passed FDA 510(k) pre-marketing notification and meets DOT/NHTSA model specifications for breath alcohol screeners.

Features

-  **Simple and easy to use**
-  **Light weight and compact**
-  **Folding mouthpiece design**
-  **Three replaceable mouthpieces (included)**
-  **Two-digit LED display**
-  **Beeps when ready to use**
-  **Beeps after breath sample obtained**
-  **Semi-conductor oxide sensor**
-  **Electric pressure sensor to improve breath capture quality**

Specifications

Dimensions	5 x 1.75 x .75 inches	Digital Display	2 Digits (.00% B.A.C.)
Weight	.3 lbs	Operating Temp	10 - 40 c (50 - 104 f)
Housing	Shock Resistant, Molded Plastic	Detection Range	.00% - .40% B.A.C.
Power	2 AA Batteries	Air Sample	5 Seconds
Battery Life	150 Tests¹	Calibration	Required Periodically³
Sensor	Sensitive Semi-conductor Sensor	Operation	1 Button
Blowing Time	5 Seconds	Limited Warranty	1 Year
Response Time	5 Seconds¹	Accuracy	+/- 0.01% at 0.02% B.A.C.²

¹ Estimates, results may vary.

² Read owner's manual before use and follow instructions carefully.

³ Visit www.q3i.com

© 2010 Q3 Innovations. All Rights Reserved.



Q3 Innovations www.q3i.com
 2349 Jamestown Avenue, Suite #4, Independence, Iowa 50644
 Tel: 319-334-3412 Fax: 319-334-3421 Email: sales@q3i.com